

## NE1 Series Unstirred Thermostatic Water Baths



The **NE1 series**, Thermostatic Clifton analogue baths are designed for general use and are economically priced without compromising quality.

Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance.

Heaters and temperature sensors are concealed under the tank providing a clutter free working area which is easy to wipe clean.

NE1-2.5		NE1-4		NE1-8		NE1-14			
GENERAL									
Temp. Range		Ambient + 5°C - 100°C							
Temp. Control		Analogue Dial Pointer							
Sensitivity/Uniformity		±0.5°C / ±0.5°C		±0.25°C / ±0.1°C					
Heater Power (220-240V)		250W		400W		800W		1000W	
Voltage		230V 50-60Hz		120V or 230V models available 50-60Hz					
DIMENSIONS									
Tank Capacity		2.5 Litre		4 Litre		8 Litre		14 Litre	
External (mm)		270w x 190d x 240h		332w x 185d x 290h		332w x 270d x 290h		361w x 332d x 290h	
Internal (mm)		230w x 130d x 150h		300w x 150d x 150h		300w x 240d x 150h		325w x 300d x 150h	
Water Level Max. (mm)		130							
Water Level Min. (mm)		55							

NE1-22		NE1-28		NE1-56	
GENERAL					
Temp. Range	Ambient + 5°C - 100°C				
Temp. Control	Analogue Dial Pointer				
Sensitivity/Uniformity	±0.25°C / ±0.1°C				
Heater Power (220-240V)	1500W			2500W	
Drain (12mm Bore)	Drain Outlet			Drain Tap Fitted	
Voltage	120V or 230V models available 50-60Hz				
DIMENSIONS					
Tank Capacity	22 Litre	28 Litre		56 Litre	
External (mm)	538w x 332d x 290h	538w x 332d x 290h		660w x 540d x 290h	
Internal (mm)	500w x 300d x 150h	500w x 300d x 200h		620w x 500d x 200h	
Water Level Max. (mm)	130	180			
Water Level Min. (mm)	55				

Accessories on Page 14-17