

NE1D Series Unstirred Digital Water Baths



NE1D-14



Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. The NE1D series, Clifton Digital Baths are designed for immersion of flasks/bottles/racks in a highly stable temperature environment throughout the water bath.

	NE1D-4	NE1D-14	NE1D-22
GENERAL			
Temp. Range	Ambient + 5°C - 99°C		
Temp. Control	Digital, LED display showing actual or set temperature, resolution 1°C		
Sensitivity/Uniformity	±0.2°C / ±0.1°C		
Heater Power (220-240V)	400W	1000W	1500W
Voltage	120V or 230V models available 50-60Hz		
Drain (12mm bore)			Drain Outlet
DIMENSIONS			
Tank Capacity	4 Litre	14 Litre	22 Litre
External (mm)	332w x 185d x 290h	361w x 332d x 290h	538w x 332d x 290h
Internal (mm)	300w x 150d x 150h	325w x 300d x 150h	500w x 300d x 150h
Water Level Max. (mm)	130		
Water Level Min. (mm)	55		

Accessories on Page 14-17